



Programação dos Reagentes CEPA

COBAS MIRA / MIRA PLUS

Uréia UV Cepa

Urease/GluDH -Cinético UV

GENERAL		CALCULATION	
MEASUREMENT MODE	ABSORB	SAMPLE LIMIT	NO
REACTION MODE	R-S	REACTION DIRECTION	DECREASE
CALIBRATION MODE	CALIBRATOR	CHECK	ON
REAGENT BLANK	REAG/DIL	CONVERSION FACTOR	1
CLEANER	NO	OFF SET	0
WAVELENGTH	340	TEST RANGE LOW	0
DECIMAL POSITION	0 mg/dL	HIGH	300
UNIT	12 mg/dL	NORMAL RANGE LOW	10
		HIGH	50
ANALYSIS		NUMBER OF STEPS	1
POST DIL FACTOR	3	CALCULATION STEP A	KINSEARCH
CONC. FACTOR	NO	READING FIRST	5
SAMPLE CYCLE	1	LAST	10
VOLUME	2 µL	REACTION LIMIT	NO
DILUTION NAME	H2O	CALIBRATION	
VOLUME	10 µL	CALIBRATION INTERVAL	ON REQUEST
REAGENT CYCLE	1	TIME	NO
VOLUME	270 µL	BLANK RANGE LOW	NO
		HIGH	NO
STANDARD POS	*	REAGENTE RANGE LOW	NO
STD 1	*	REAGENT RANGE HIGH	NO

*Inserido pelo usuário

Para Cobas Mira Plus, substituir CALIBRATOR por SLOPE AVG