

Creatinina Cepa

Jaffé-Cinético de Tempo fixo

TEST TYPE	Fixed Time Single
SERUM STARTER	INACTIVE
FILTERS A/B	510/700
UNITS	mg/dL
TEST METHOD	WHIT FACTOR
TEST METHODOLOGY	
NUMBER OF WASHES	2
DELAY TIME	0
INC REAGENT	30
READING TIME	60
TEST LIMIT (CONC)	20
MAX ABS DELTA	200
REAGENT mABS Limit	600
REAGENT A/B	300
REACTION DIRECTION	Increasing
REAGENT DILUTION	01:01:00
INITIAL ABS	1800
CURVE ACCEPTANCE	100,00%
AUTOMATIC PROFILE	(INACTIVE)
RERUN TEST RGT BLK H:M	00:00:00
SERUM PARAMETERS	
TEST'S NAMES	Creatinina
SAMPLE VOLUME	30
DILUTION RATIO	1:3
MIN. MAX. M	0.6 – 1.1
MIN. MAX. F	0.5 – 0.9
MIN. MAX. B	