



Programação dos Reagentes CEPA

RA-XT 500 - 1000

Triglicérides Cepa

Enzimático Monoreagente

CHEM #	(*)
NAME	TRIGLICÉRIDES
IMMUNOASSAY	NO
TYPE	ENDPOINT
INVERSE	NO
% SAMPLE	6 (3µL)
FILTER P	500
BIC. CHEM	NO
DELAY TIME	4
BLANK	NO
%REAG. VOL	60 (300µL)
2° REAG.	NO
UNITS	mg/dL
UNIT FACTOR	1
DECIMAL	0
RBL LOW	0
RBL RIGH	0,05
CAL FACTOR (**)	*
STANDARD VALUE	*
NORMAL LOW	60
NORMAL RIGH	150
SLOPE	1
INTERCEPTION	0
END POINT LIMIT	0,02
C1 *	*
D1 *	*
DELTA#	*