



Programação dos Reagentes CEPA

Humastar 80

ASLO TURB CEPA

Turbidimetria

Description	Anti-streptolisina O	
Test ID	ASLO TURB CEPA	
Manufactured	CEPA	
Expire	-	
Modo	End Point	
Note	-	
Wavelength (nm)	540	
Volumes (µL)	Sample	5
	Reagent 1	350
	Reagent 2	-
Reading parameter	1 st incubation	10
	2 nd incubation	120
	Stability	0
	No Reading	2
Measure units	UI/mL	
Decimal Point	1	
Minimum conc.	1	
Calibration	Specific	
N Standard	1	
Replicate	1	
Offset	0	
Concentration	*	
Function	Incr.	

*Vide Bula